LIPOCUBE

TR ÜRÜN KULLANIM KILAVUZU

EN INSTRUCTION FOR USE

TR

Lipocube PRP Kiti 1
Sodyum Sitratlı Steril PRP (Plateletten Zengin Plazma) Kan Seperasyon Kitleri
Referans no: OM-CU1-0001

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PRODUCT DESCRIPTION and PURPOSE

With Blood Separation Kits, the patient's blood components are separated according to their density using centrifugation. Thanks to the better quality PRP can be obtained from the blood coagulation of the sample taken.

PACKAGE CONTENT

2 pcs 10ml Sodium Citrate (1ml) Vacuum tube
1 pc 5 ml luer lock injector
1 pc 21G 0.8*40 mm Needle
1 pc 18G 90 mm Needle
1 pc 21G 0.8*38 mm Needle
1 pc Holder

INDICATIONS

The product is indicated in the treatment of osteoarthritis and lateral epicondylitis. In addition, PRP can be used as a supplement to standard treatments in cases of atopic dermatitis, alleviation of disease symptoms and chronic wounds, lupus (butterfly disease), alopecia areata / barbae / androgenetic diseases.

CONDITIONS

PRP application is not recommended in the presence of any of the following conditions:

- Platelet Dysfunction Syndrome
- Critical Thrombocytopenia
- H, pofibrinogenaemia
- Hemodynamic instability
- Autoimmune diseases
- Malignancy
- Sepsis
- Acute chronic infections
- Chronic Liver Pathology
- Anti-coagulation therapy (ASA therapy, etc.)
- Pregnancy
- Newborns

WARNINGS and PRECAUTIONS

1. The kits are unique to each patient and are for single use only.
2. It must be ensured that the additional materials (injector, needle tip) to be used are reliable and appropriately certified.

Document No: TD.01.121-8 | Release Date: 10.11.2021 | Revision No:00 | Revision Date: -

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2 years from the date of manufacture.

USAGE

- Open the box.
- Open the sterile blister pack.
- Combine the holder and the blood collection needle.
- Apply the blood collection needle to the cleared area where blood will be taken from the patient's arm. (Venous)
- Place the PRP tube into the holder where the needle is attached. (Routine Phlebotomy)
- PRP Tube is vacuumed. Once the tube receives sufficient blood, it stops on its own. (Do the same procedure for both tubes without removing the blood collection needle from the vein.)
- Centrifuge the tube in the centrifuge device at 2500 - 4000 RPM for 8 minutes. Centrifuge gots must be compatible with tubes.
- Attach the 18G 90 mm needle to the 5ml Luer-lock syringe.
- Collect PRP over the red blood (RBC) layer by inserting an 18G 90 mm needle into the PRP tube. (Do not open the cap of the tubes unless Laminar Flow Conditions are met.)
- Collect the 2ml PRP shot in a 5ml syringe with caution.
- There will be 2ml of PRP in the syringe.
- Gently shake the syringe for 60 seconds to homogenize possible cell clumps in the plasma.
- Attach the needle to be applied to the 5ml Luer-Lock syringe.
- Apply to the patient by injection method.

Explanation of Symbols on Box Label

No	Symbol	Description
1		Producer
2		Expiration Date
3		Refer to the user manual.
4		Caution, read the user manual before use.
5		Do Not Reuse
6		LOT Number
7		Radiation Sterilization

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8		CE Symbol of SZUTEST Test Notified Body
9		Do not use damaged or opened packaged products.
10		Temperature Limits
11		Keep dry.
12		Keep away from sunlight.
13		Do not re-sterilize.
14		Model Name
15		Production Date
16		Biohazard
17		Ethylene Oxide Sterilization

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Endüstri ve Sanayi Dış Ticaret Limited Şirketi