TR

Lipocube PRP Tüpleri  
Sodyum Sitratlı Steril PRP (Plateletten  
Zengin Plazma) Kan Seperasyon  
Tüpleri Referans no: OM-L1-0006

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## PRODUCT DESCRIPTION and

### PURPOSE

With Blood Separation Kits, the patient's blood components are separated according to their density using centrifugation. Thanks to the an better quality PRP can be obtained from the blood coagulation of the sample taken.

### PACKAGE CONTENT

6 pcs 10ml Sodium Citrate (1ml) Vacuum tube Holder

### INDICATIONS

The product is indicated in the treatment of osteoarthritis and lateral epicondylitis. In addition, PRP can be used as a supplement to standard treatments in cases of atopik dermatitis, alleviation of disease symptoms and chronic wounds, lupus (butterfly disease), alopecia areata / barbae / androgenetic diseases.

### CONDITIONS

PRP application is not recommended in the presence of any of the following conditions:

- Platelet Dysfunction Syndrome
- Critical Thrombocytopenia
- H<sub>2</sub> pofibrinogenaemia
- Hemodynamic instability
- Autoimmune diseases
- Malignancy
- Sepsis
- Acute chronic infections
- Chronic Liver Pathology
- Anti-coagulation Therapy (ASA therapy, etc.)
- Pregnancy
- Newborns

### WARNINGS and PRECAUTIONS

1. The kits are unique to each patient and are for single use only.  
2. It must be ensured that the additional materials (injector, needle tip) to be used are reliable and appropriately certified.

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3. For medical applications only. It is not used for other purposes.
4. It should only be used by medically licensed physicians.
5. Expired products should not be used.
6. Do not use products that have been opened or damaged.
7. Calibration of the centrifuge device to be used must be done. If a fixed-angle centrifuge must be used, a centrifuge with a minimum opening of 45 degrees must be used.
8. PPP should be used immediately after PRP has been separated.
9. The patient should be informed about the risks and possible side effects of the treatment.
10. If anti-inflammatory drugs and blood thinners are used, they need to be discontinued at least one week before the treatment. The use of anti-inflammatory drugs and blood thinners is not recommended for the 2 weeks following the treatment.
11. During the day of treatment, excluding hot bath, sauna, and doctor's recommendations, no heavy exercise should be performed.
12. Consideration should be given to the amount of samples in case of vacuum loss.

## PATIENT POPULATION

The product should be used on adult patients.  
User Profile:

The individuals who will use our device must;  
a) Be professionals who have completed all their medical education,  
b) Be careful during the operation,  
c) Have sufficient experience in the use of the device.

## SIDE EFFECTS

Pain in the application area, temporary bruising, or edema may occur in skin applications. Vascular damage, hematoma, and/or infection may occur if the product is applied incorrectly by non-experts and sterility conditions are not followed.

## STERILIZATION

PRP tubes were sterilized by Gamma Sterilization and set Ethylene Oxide method. It is disposable. It cannot be re-sterilized and used.

## STORAGE CONDITIONS

Must be kept in a special compartment in a clean environment at a temperature of 15-25 degrees. Keep away from damp and direct sunlight areas. Do not place heavy objects on it.

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



## SHELF LIFE

2 years from the date of manufacture.

No	Symbol	Description
1		Producer
2		Expiration Date
3		Refer to the user manual.
4		Caution, read the user manual before use.
5		Do Not Reuse
6		LOT Number
7		Radiation Sterilization
8		CE Symbol of SZUTEST Test Notified Body
9		Do not use damaged or opened packaged products.
10		Temperature Limits
11		Keep dry.
12		Keep away from sunlight.
13		Do not re-sterilize.
14		Model Name
15		Production Date
16		Biohazard
17		Ethylene Oxide Sterilization

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Endüstri ve Sanayi Dış Ticaret Limited Şirketi